

THE *SMART WAY* TO GET *THINGS DONE*




Cappello™ Automated Vial Cap Dispenser

***CONTINUOUS FEED – COMPACT – QUIET
FULLY ELECTRIC (no compressed air required)
USES PRE-LOADED CAP CASSETTES***

***The Cappello™ automated vial cap dispenser
can be integrated with any automated laboratory
workstation where fresh caps are needed.***

- Small, Compact and Quiet (replaces large noisy bow feeders).
- Pre-loaded cap cassettes save time and money.
- Easy, fast reload.
- Cap cartridge housing holds multiple cassettes for large cap capacity before reload required (>500 cap capacity for cryovials and glass vials).
- Continuous feed for fast, fully automatic dispensing.
- Fully electric operation eliminates need for compressed air.
- Easy to operate.

Dimensions: 5"h x 11"w x 12"l



Innovative Laboratory Products 

1.888.323.3585

www.tradewindsdirect.com

LAKEVIEW CORPORATE PARK • 10555 86TH AVENUE • PLEASANT PRAIRIE, WI 53158 • FAX: 262.605.3262