

THE *SMART WAY* TO GET *THINGS DONE*

High Recovery Vials & Tubes

TWD high recovery vials and tubes are designed for maximum powder and liquid transfer recovery

- Minimal Sample Loss
- Easy Transfers
- Standard Exterior Dimensions



96 well glass high recovery micro tubes with ceramic/laser-etched 2D barcodes

Glass High Recovery 96-Well Microtubes

Ceramic/Laser-etched 2D and/or 1D Barcodes

Flanged opening is designed to use standard caps and fit into commercially available 96-well racks.

- Optional laser etched 2D barcodes on microtube bottom and/or linear bar codes on tube side.
- Black on white or white on black
- 2 standard heights
- Custom heights available



GHRMT-7.5X42

7.5mm x
41.5mm -
1.0mL

GHRMT-7.5X48

7.5mm x
48mm -
1.2mL



High Recovery

LIQUID TRANSFER VIALS

Steeply sloped inner conical surface to ensure maximum recovery of liquid samples. Same outer dimensions and threads as standard 1-dram, 2-dram, and 4-dram vials.



Innovative Laboratory Products

1.888.323.3585

262.605.3260 • www.twdinc.us

LAKEVIEW CORPORATE PARK • 10555 86TH AVENUE • PLEASANT PRAIRIE, WI 53158 • FAX: 262.605.3262

THE *SMART WAY* TO GET *THINGS DONE*

High Recovery Vials and Tubes



High Recovery



SOLID/LIQUID TRANSFER VIALS

Smoothly sloped inner conical surface ensures a high level of solid/liquid recovery. Dead volume <50ul as measured using a 20G hyperdermic needle.

STANDARD SIZES

- HR1232C** - 12x32mm 1mL/autosampler vials
 - HR1545C** - 15x45mm, 4mL/1 dram, 13-425 caps
 - HR1760C** - 17x60mm, 8mL/2 dram, 15-425 caps
 - HR2857C** - 28x57mm, 20mL/4 dram, 24-400 caps
 - HR2895C** - 28x95mm, 40mL/8 dram, 24-400 caps
- Custom sizes, clear/amber glass



Straight-Wall High Recovery

SOLID/LIQUID TRANSFER VIALS

Smoothly sloped inner conical surface ensures a superior level of solid/liquid recovery. Dead volume <10ul as measured using a 20G hyperdermic needle.

STANDARD SIZES

- SWHR12.7X45** - 12.7x45mm, 13-425 caps
 - SWHR14X60** - 14x60mm, 15-425 caps
- Custom sizes, clear/amber glass



TWD TradeWinds, INC

Innovative Laboratory Products

1.888.323.3585

262.605.3260 • www.twdinc.us

LAKEVIEW CORPORATE PARK • 10555 86TH AVENUE • PLEASANT PRAIRIE, WI 53158 • FAX: 262.605.3262