

## SPE-Columns

### with Pierceable Cap

Polypropylene SPE column with Open-top screw cap & 5-10micron PE frits.  
Frit material can be made to specifications.

#### Description

- Polypropylene SPE column:
- Dimensions with caps - 73mmh x 35.5mm O.D
- Maximum capacity – 25mL
- Frits –5-10micron PE frit disk
- Cap size: 34-400 GPI, open top with sil/TEF septa.
- Caps are packed separately



150 mm long • Polypropylene

### DISPOSABLE

## Polypropylene Test/Culture Tubes

### Extra-Long • Tall

Part #	Polypropylene Test/ Culture Tube Size (mm)	Maximum Capacity (mL)	Items/pk
DPT16x150	16x150	15	250
DPT18x150	18x150	25	250
DPT20x150	20x150	35	250
DPT25x150	25x150	55	250

## Glass Round-Bottom Vials

- The round bottom vials are designed to improve transfer recovery.
- All vials are made of clear USP type I borosilicate glass and use standard 24-400 caps.

Part #	Size (mm)	Capacity (mL)	Items/pk	# pk/cs
GRV-20	28x57	20	120	6
GRV-20TALL	28x70	27	144	5
GRV-40	28x95	40	144	3
GRV-40TALL	28x108	47	120	3

Contact your TWD Representative for Information On

## Hi-Recovery Vials



Increase  
Transfer  
Recovery



Innovative Laboratory Products 

**1.888.323.3585**

262.605.3260 • [www.twdinc.us](http://www.twdinc.us)

LAKEVIEW CORPORATE PARK • 10555 86TH AVENUE • PLEASANT PRAIRIE, WI 53158 • FAX: 262.605.3262